

MCP IMAGE INTENSIFIERS



FEATURES


- Custom design
- Proximity focused - No distortion
- 12, 18, 25, 40 and 75mm
- High gain options
- Compact size
- X-ray, VUV, solar blind, visible and IR
- Fast gating
- Unaffected by magnetic fields

APPLICATIONS INCLUDE

- Spectroscopy
- Low light level imaging
- Astronomy
- Photon counting
- Multi-channel spectrophotometry
- High speed imaging
- Biofluorescence
- Bioluminescence
- Time resolved imaging

INTRODUCTION

Photek manufactures a range of 12, 18, 25, 40 and 75mm active diameter image intensifiers. The metal ceramic body is rugged and the proximity focus design gives a distortion free image in a very short overall length. A range of photocathodes and input window materials enables a wide choice of spectral responses to suit many applications. The fibre optic output ensures a defined output focal plane and allows efficient coupling to CCDs and linear image sensors. A variety of MCP configurations satisfies all gain requirements.



Spectral Response and Phosphor Screen

Photek are able to offer a wide range of Photocathodes including Solar Blind, Bialkali, LNS20 and S20. Phosphors include P11, P20, P24, P31, P43 P46 and P47.

Please refer to the separate Photocathode and Phosphor screen data sheets for further details.


Input Windows

All Photek intensifiers are available with a choice of input window materials. These include Glass, Fibre Optic, Fused Silica and MgF₂.

Resolution

The resolution of an image intensifier is reduced as the number of MCPs is increased. The modulation transfer function (MTF) measures the relative output contrast as a function of spatial frequency when a suitable bar pattern is imaged by the tube under test.

Diode image converters without an MCP have almost no gain, but can achieve high MTF and a limiting resolution of 50 lp/mm. They are useful as optical shutters and wavelength converters.



The superior pulse height distribution is used to discriminate between photon events and electronic noise, and the resolution performance can be recovered by suitable event-processing software. Photon counting systems can in fact achieve resolution only limited by the size of the microchannel plate pores (usually 10 microns).

Uniformity

The standard overall uniformity of Photek intensifiers is:

Detector Size	Uniformity
12mm	To be determined
18mm	+/-15% (Typ +/-10%)
25mm	+/-20% (Typ +/-12%)
40mm	+/-25% (Typ +/-15%)
75mm	+/-30% (Typ +/-17%)

Photographic Applications


In photographic applications it is common to set the tube up in contact with a film and switch all voltages, except the photocathode, to normal running conditions. It is important that no light is produced by the tube until the photocathode is switched on to make the photographic exposure. Our film grade tubes are tested for exposures of up to 10 min in contact with ASA 3000 film.

Space Qualified and Other Special Image Intensifiers

Photek is involved with several space projects, and has experience in both the mechanical and quality assurance aspects of these programmes. We have developed special ruggedised structures, housings and space encapsulation techniques for these projects.

Gain


The MCP voltage is used to control the gain, and each tube is supplied with test data showing gain as a function of MCP voltage. Photek routinely manufactures a range of image intensifiers with one, two or three MCPs (both single and double thickness), and with a variety of configurations. Photek also manufactures proximity photodiodes using the same construction but without any MCPs.



Pulse Height Distribution

Single MCP intensifiers have an exponential pulse height distribution.

For photon counting image intensifiers, pulse height distribution is measured during processing to optimise channel plate outgassing.



Gating/Electronic Shuttering

Many customers use these tubes as fast electronic shutters. Applying a voltage pulse to the photocathode turns a proximity focus intensifier into an optical block. The on/off gain/attenuation ratio is typically greater than 10¹⁰. Tubes can be supplied with a small gating control unit which will accept a 5 Volt TTL pulse to turn the tube on/off in less than a microsecond, and rather more sophisticated units can achieve control down to a few ns.

For fast gating response it is usually necessary to increase the photocathode conductivity by using a transparent conducting undercoat or by reducing cathode size. Photek make tubes of both types, the most common of the latter type being the slot photocathode for multi channel analyser applications. 25mm intensifiers can be switched on/off in approximately 3ns FWHM using these techniques, while 40mm tubes have achieved 12ns FWHM. Tubes can also be modulated at high frequency; over 100MHz has been achieved with a special 25mm intensifier.

Ultra-Fast Gating

Ultra-fast gating is possible using a conductive mesh undercoat. The inherent iris time can be reduced to 50ps by this technique, but there are no commercially available pulse generators made to achieve this. Electrical connection to the tube is also a problem; Photek has experience of building tubes with stripline and other connector designs which have

Blemishes

The standard Photek blemish specification is given below.

Zone (mm)	18	18 to 25	25 to 40	40 to 75
75 to 100	3	5	9	50
101 to 150	2	3	6	20
151 to 250	0	2	3	8

Photek Ltd

26 Castleham Road,
St Leonards on Sea,
East Sussex, TN38 9NS,
United Kingdom

Tel (+44) 1424 850555
Fax (+44) 1424 850051
Email sales@photek.co.uk
Website http://www.photek.co.uk

Photek Inc.,

2030 East Broadway
Suite 206,
Tucson, AZ 85719,
USA

Tel (520) 622 1289
Fax (520) 622 0893
Email photekinc@compuserve.com
Website http://www.photek.com